



Secura®

<b>Secura 26P-10IN Hallmarking Balance</b>
--

**Specifications**

Model		SECURA26P-10IN
Accuracy class		<b>Ⓛ</b>
Type		SQP-H
Max. capacity	g	1.1   2.1   3.1   21
Scale interval d	mg	0.001   0.002   0.005   0.01
Verification scale interval e	mg	1
Min.	mg	0.1
Repeatability (standard deviation)	≤ ±mg	0.003   0.004   0.005   0.03
Repeatability (standard deviation), typical	≤ ±mg	0.0025   0.003   0.004   0.02
Linearity deviation	≤ ±mg	0.01
Sensitivity drift between +10 to +30 °C	± ppm/K	1
Tare equalization range (subtractive)		< 100% from max. weighing capacity
Typical stabilization time	≤ s	6   6   5   5
isoCAL: – Temperature change – Time interval	K h	1.5 4
Display result (depending on the set filter level)	s	0.2 - 0.4
Weighing pan size	mm	∅ 50
Weighing chamber height*	mm	218
Net weight, approx.	kg	8.0

\* Upper edge of the weighing pan to the lower edge of the upper draft shield panel